



ELEVAND STEADYSTAND



ELEVAND STEADYSTAND

The Elevand SteadyStand enables a single caregiver to assist patients or residents with sit-to-stand transfers throughout the day. This device promotes early mobilization by encouraging patients to be active during these transitions. Designed to accommodate a range of patients with different clinical needs, the Elevand SteadyStand offers versatility in a single size. Additionally, it is non-powered, providing a simplified and efficient solution for everyday transfers.

FEATURES

Seat Parts

Rotate sideways, making it simpler for the end-user to get into position.

Knee Pad

Provides essential support for raising the user into a half-standing position and stabilises the patient's legs during the transfer.

Rear Castors with Lock:

Prevent unintentional movement, ensuring stability while the end-user is sitting down.



Two Large Handles

One for the assistant to manoeuvre the device and another for the end-user to hold onto when entering the lift or during transfer.

Foot Lever

Allows for easy adjustment of the base leg spread, and when combined with the low footrest, facilitates easy access to wheelchairs, beds, and other surfaces.



Maximum Safe Working Load
180kgs

TECHNICAL SPECIFICATIONS

Base leg spreading	Mechanical
Max. safe working load (SWL) (Kg)	180
Length (cm)	90.5
Width (cm)	62
Height (cm)	102
Knee pad height (cm)	64
Min.base legs width (int./ext.)(cm)	49.3/61.1
Max.base legs width (int./ext.)(cm)	91.8/103.6
Base legs height (cm)	8.4
Castors diameter (without/with lock) (cm)	7.5/10
Weight (kg)	35.5
Min.turning diameter (cm)	110
Packed dimensions (LxWxH)(cm)	120 x 80 x 110

AVAILABLE DOWNLOADS



 [Instruction for Use](#)

 [Quick Start Guide](#)

 [Technical Specifications](#)